

User Guide


Wi-Fi Smart Automatic Pill dispenser

M228A



Product Overview


1.Product overview




Manual*1




Key*2



Charge cable*1




Adapter*1



Dosage ring*9

Product Overview

2.LCD Display Overview



** Note: The product, accessories, user interface, and other illustrations in this quick start guide are for schematic purposes only and are for reference. There may be slight differences between different versions. Please refer to the actual product.*

Product Overview

3.Touch Button

Icon	Name	Quick Press	Long Press
	Take Pills Button	Take Pills	Early Dose
	Manual Button	Enter Manual /Determine Options	RestartCycle
	Back Button	Back	Child Lock ON/OFF
	SOS Button		This version is not supported.
	Roll Up Botton	Roll Up	
	Roll Down Button	Roll Down	

**Note: If the touch panel is not operated for a period of time, it will automatically lock. Please remember to long-press the  key to unlock.*

**Note: The [Emergency Call] function is supported in some versions. Please refer to the actual product. This function can only be used when the device is connected to the network and requires the payment of a certain service fee.*

Operating Guide

4. Power Supply


Please plug in the adapter to the device and open the power supply switch

***Inner battery is only for back up, please keep in charging during the daily use.**

5. Device binding

Smart pill box devices need to be connected to the Internet to bring users a better experience. Please use your smartphone to configure the network for the smart pill box.

- 5.1 Scan the QR code to download and install the app on the Phone.
- 5.2 Open the App, log in/sign up the account.
- 5.3 Binding: In the [Home] of the App, click [Add], scan the QR code at the bottom of the device.
- 5.4 Edit the device name and family member to bind the device.



Operating Guide

5.5 Unbind the device, you can click the [Settings] button in the upper - right corner of the APP, and select [Authorization and Unbinding] - [Unbind].

5.6 Exit the WiFi connect operating: After the device finishes network configuration, it will automatically exit the network - configuration mode. If you need to force an exit, please press the ↵ , and the device will be forced to exit the network configuration mode.

6.Setting Medication Schedule

After loading the medicines, setting the "Active Slots", and "Alarms". The APP will be updated and upgraded in the future, and the interface may be changed, but it will not affect the user's use.

Operating Guide


7. Load the Pills

7.1 Use the provided key to unlock the pillbox, and press the lock core to lock the pillbox.




7.2 Lightly press the pill tray and rotate it counter-clockwise to remove the pill tray.

7.3 Rotate the pill tray lid counter-clockwise to remove it.



7.4 Please rotate the pill tray till the tray stuck.

7.5 According to the doses taken per day, configure the dosage rings and place the dosage ring's positioning piece in the starting slot



Operating Guide

7.6 Please follow the instructions on the dosage rings and load the pills counterclockwise.

7.7 After loading the pill tray into the main structure, you need to proceed the "Restartcycle" operation

***Note: Refer to standalone operation 10.8 or use the "RestartCycle" function in the app.**


8. Take Pills

While the alarm rings, the user needs to press  to take the pills. After the pills falls into the medicine cup, the user takes out the medicine cup and takes the pills.

Users can also use the "Take in advance" function to take pills of next alarm. Press and hold  to rotate the pill tray in advance. Meanwhile, press  again to take the pills.

a. Press 

b. Take out the medicine cup



Operating Guide

9. Operating the Dispenser

Press  to enter the manual, press  and  to roll, choose the options, press  to choose, press  to cancel or back to last level.

9.1 Time Setting

Select "Clock" on the Manual, press  and  to adjust the Year, Month, Day, Hour, Minutes, Press  to confirm each selection.

Wi-Fi connected status will auto sync time.

9.2 Alarms Setting

Select "Alarm" on the Manual, press  and  to scroll up and down respectively between the alarm slots. Then select a slot with , press  and  to adjust hour, press  to confirm, press  and  to adjust minutes, press  to confirm the selection, setting all alarms by the same way. Press  to exist alarm setting

Cancel alarms, select the alarm, long press  to clean the alarm.

9.3 Date Format, Time Format, Volume, Language, Take advance, Ring Mode, Ring Type, Dispensing setting

Use  and  to select the options from the manual, press  to confirm.

**Regarding Take advance:

Disable: It is not possible to dispense medicine in advance on the device.

Enable: Long press  to dispense medicine in advance.

Lock: The child lock needs to be unlocked before the medicine can be dispensed in advance.

Double Check: After long press , you need to press  again to dispense the medicine in advance.

**Regarding Dispensing:

Button: When the medication time arrives, users need to press  for the medicine to fall into the medicine cup.

Auto: When the medication time arrives, the medicine will automatically fall into the medicine cup.

9.4 Alarm effective period

Operating Guide

Turn on this function, you can set the alarm effective period to manage the med course for the user.

Select "Med Course" from the manual, press ▲ and ▼ to select numbers of year, month, day of the start time, press 📅 to confirm your selections, same way to confirm the end time.

Long press 📅 to clean the med course period at same interface.

9.5 Mute Set

During time period of silence mode, the device is in mute status, when it's time for medication, the device will remind the user by the flash light.

Select "Mute Set" on the manual, press 📅 to switch between "Mute On" and "Mute Off"

Use ▲ and ▼ to adjust hour and minute, press 📅 to confirm the start time and end time.

9.6 Sound Record

The user can record the sounds by themselves to replace the original reminding voices on the device.

1. Select "Sound Record" on the manual, use ▲ and ▼ to scroll up and down.
2. Long press 📅, after hearing a beep sound, start the sound record, during time of recording, please hold on the 📅 button till finish the recording operation.
3. Move away the fingers from the 📅 button after finishing the record, original voices has been replaced to new sounds.
4. Press 📅 once to broadcast the sound to confirm the effect.

***This device support for multiple times record till the user is satisfied.**

"Reset All"-Clean all sounds recorded by the users, reset to original reminding sounds.

9.8 RestartCycle

After loading the pill tray, please don't forget to "RestartCycle", use ▲ and ▼ to scroll up and down to select "RestartCycle" on the manual, press 📅 to confirm, select "Confirm".

9.9 Factory Test

Use ▲ and ▼ to select the options from the manual, press 📅 to confirm. Enter this interface, test the functions.

***It's special for factory test.**

Trouble Shooting

Trouble shooting

Icon	Problem and Possible Solutions
	The device is abnormal: Please contact the dealer or factory to processing.
	Pill tray rotate failure: Please check if there are any foreign stuffs stucked on the tray, remove the stuffs, and try to restart cycle.
	Restart cycle failure: Please check if there are any foreign stuffs stucked on the tray, remove the stuffs, and try to restart cycle.
	Touch button failure: Please contact the dealer or factory.
	Door of the pill tray failure: Please check if there are any foreign stuffs stucked the door, load the pill tray again after removing the stuff.
	The push rod of door failure: Please contact the dealer or factory to processing.
	Refilling remind: The pill tray is empty, please refill the pills in time, or check if the active slot number is correct or not on the app.
	T&H Alarm Temperature and humidity alarm is displayed when the internal temperature and humidity of the device exceeds the set value (temperature: <0°C, >38°C humidity: >70%). Please pay attention to the environment where the device is placed.